



# Variant: Released

## Integrated Actuators Laboratory (LAI)

Semester project: 02/2025 - 06/2025

Student: Théo Heng - theoheng@icloud.com

Supervisors: Maël Dagon, Paolo Germano

2025-05-27

Rev 1.0

| Page | Index                           | Page | Index                    | Page | Index                             | Page | Index            |
|------|---------------------------------|------|--------------------------|------|-----------------------------------|------|------------------|
| 1    | Cover Page                      | 11   | Motor Control - Drivers  | 21   | Motor Control - Driver 9          | 31   | Revision History |
| 2    | Block Diagram                   | 12   | Motor Control - Driver 0 | 22   | Motor Control - Driver 10         | 32   | .....            |
| 3    | Project Architecture            | 13   | Motor Control - Driver 1 | 23   | Motor Control - Driver 11         | 33   | .....            |
| 4    | MCU - Unit                      | 14   | Motor Control - Driver 2 | 24   | Motor Control - Driver 12         | 34   | .....            |
| 5    | Power - Battery Management      | 15   | Motor Control - Driver 3 | 25   | Motor Control - Driver 13         | 35   | .....            |
| 6    | Power - Wireless Receiver & USB | 16   | Motor Control - Driver 4 | 26   | Motor Control - Driver 14         | 36   | .....            |
| 7    | Power - Path & Generation       | 17   | Motor Control - Driver 5 | 27   | Motor Control - Driver 15         | 37   | .....            |
| 8    | Logic - I2C Mux                 | 18   | Motor Control - Driver 6 | 28   | Motor Control - Driver 16         | 38   | .....            |
| 9    | User - RGB LED Indicator        | 19   | Motor Control - Driver 7 | 29   | Motor Control - Driver 17         | 39   | .....            |
| 10   | Sensing - IMU                   | 20   | Motor Control - Driver 8 | 30   | Misc - Mounting Holes & Fiducials | 40   | .....            |

## TOP VIEW



## BOTTOM VIEW



## DESIGN CONSIDERATIONS

DESIGN NOTE:  
Example text for informational design notes.

DEBUG NOTE:  
Example text for debug notes.

DESIGN NOTE:  
Example text for cautionary design notes.

DESIGN NOTE:  
Example text for critical design notes.

LAYOUT NOTE:  
Example text for critical layout guidelines.

## NOTES

### MULTIPLE MOTOR CONTROLLER WITH SENSORLESS POSITION CONTROL

Not fitted components are marked as

DRAFT - Very early stage of schematic, ignore details.

PRELIMINARY - Close to final schematic.


CHECKED - There shouldn't be any mistakes. Contact the engineer if you find any.

RELEASED - A board with this schematic has been sent to production.

Released - 24/04/2025

|   |  |                        |                                       |
|---|--|------------------------|---------------------------------------|
| Comments:<br>KiCad Template by Vincent Nguyen | Laboratory:<br>Integrated Actuators Laboratory | Variant:<br>EPFL       |                                       |
| Board Name:<br><b>Plume</b>                   | Plume  | Project Name:<br>      | Smart Footwear for Diabetic Foot Care |
| Sheet Title:<br>Cover Page                    | File Name:<br>Plume.kicad_sch                  | Designer:<br>Théo Heng | Date:<br>2025-04-24                   |
| Sheet Path:<br>/                              | Supervisor:<br>Maël Dagon, Paolo Germano       | Size:<br><b>A3</b>     | Sheet:<br><b>1</b> of <b>31</b>       |

## [2] Block Diagram



| JLC081611-7628(Special/Finished thickness 1.59mm±10%) |                             | JLC081611-2116(Standard /Finished thickness 1.65mm±10%) |         |                        |             |                       |                           |
|---|-----------------------------|---|---------|------------------------|-------------|-----------------------|---------------------------|
| Impedance (Ω)   | Type                        | Signal Layer  | Top Ref | Bottom Ref             | Trace Width | Trace Spacing         | Impedance trace to copper |
| 50  | Single Ended (Non coplanar) | L1  | /       | L2                     | /           | /                     | /                         |
| <b>Layer</b>  |                             | <b>Material</b>   |         | <b>Thickness (mil)</b> |             | <b>Thickness (mm)</b> |                           |
| L1  |                             | Outer Copper Weight1oz                                  |         | 1.38                   |             | 0.0350                |                           |
| Prepreg   | 2116 RC54% 4.9mil           |   |         |                        | 4.29        |                       | 0.1090                    |
| L2  |                             | Inner Copper Weight                                     |         | 1.18                   |             | 0.0300                |                           |
| Core  | 0.25mm 1/1OZ without copper |   |         |                        | 9.84        |                       | 0.2500                    |
| L3  |                             | Inner Copper Weight                                     |         | 1.18                   |             | 0.0300                |                           |
| Prepreg   | 2116 RC54% 4.9mil           |   |         |                        | 4.29        |                       | 0.1090                    |
| L4  |                             | Inner Copper Weight                                     |         | 1.18                   |             | 0.0300                |                           |
| Prepreg   | 2116 RC54% 4.9mil           |   |         |                        | 4.29        |                       | 0.1090                    |
| L5  |                             | Inner Copper Weight                                     |         | 1.18                   |             | 0.0300                |                           |
| Prepreg   | 2116 RC54% 4.9mil           |   |         |                        | 4.29        |                       | 0.1090                    |
| L6  |                             | Inner Copper Weight                                     |         | 1.18                   |             | 0.0300                |                           |
| Core  | 0.25mm 1/1OZ without copper |   |         |                        | 9.84        |                       | 0.2500                    |
| L7  |                             | Inner Copper Weight                                     |         | 1.18                   |             | 0.0300                |                           |
| Prepreg   | 2116 RC54% 4.9mil           |   |         |                        | 4.29        |                       | 0.1090                    |
| L8  |                             | Outer Copper Weight1oz                                  |         | 1.38                   |             | 0.0350                |                           |

### Current Recommended Stackup:


8 layer, 1oz out/1oz:  
JLC081611-2116 (Finished thickness 1.65mm±10%)

### Target specifications:


|                          |               |
|--------------------------|---------------|
| Battery Input voltage:   | 3.08 - 4.25 V |
| Max charge current:      | 455 mA        |
| Nb of controlled motors: | 18            |
| Max load current:        | 4.54 A        |
| Standby Current:         | 50 mA         |
| Sleep Current:           | 15 $\mu$ A    |

|                                |   |                          |
|--------------------------------|---|--------------------------|
| Comments:                      | Laboratory:<br><b>Integrated Actuators Laboratory</b>         | Variant:<br>Released     |
| Board Name:<br><b>Plume</b>    | Project Name:<br><b>Smart Footwear for Diabetic Foot Care</b> |                          |
| Sheet Title:<br>Block Diagram  | File Name:<br>Block Diagram.kicad_sch                         | Designer:<br>Théo Heng   |
| Sheet Path:<br>/Block Diagram/ | Supervisor:<br>Maël Dagon, Paolo Germano                      | Date:<br>2025-04-24      |
|                                |   | Revision:<br>1.0         |
|                                |   | Size:<br><b>A3</b>       |
|                                |   | Sheet:<br><b>2 of 31</b> |

# [3] Project Architecture



## [4] MCU - Unit



# [5] Power - Battery Management

Current analysis:

Normal:  
BQ29737: 4uA  
MCP73831: 100nA

Shutdown:  
BQ29737: 100nA  
MCP73831: 100nA

DESIGN NOTE - gate-source resistors:

TI recommends placing a high impedance 5MO across the gate source of each external FET to deplete any charge on the gate-source capacitance.

DESIGN NOTE - EXTERNAL PROTECTION FETs:


The external FET selection is important to ensure the battery pack selection is sufficient and complies to the requirements of the system.

BQ29737 charge overcurrent threshold: -50mV  
BQ29737 discharge overcurrent threshold: 100mV

The CSD85302L when driven at 3.7V at the gate, have a source to source resistance of ~22 mΩ.

Resulting charge overcurrent protection:  
50mV / 22mΩ ≈ 2.27A

Resulting discharge overcurrent protection:  
100mV / 22mΩ ≈ 4.54A






**Voltage and Current Protection Integrated Circuit for Single-Cell Li-Ion and Li-Polymer Batteries and Miniature Single-Cell, Fully Integrated Li-Ion, Li-Polymer Charge Management Controller**

|               |       |               |
|---------------|-------|---------------|
| COUT          | TP501 | COUT          |
| DOUT          | TP502 | DOUT          |
| CHARGE_STATUS | TP503 | CHARGE_STATUS |
| BAT_SENSE_EN  | TP504 | BAT_SENSE_EN  |
| BAT_SENSE     | TP505 | BAT_SENSE     |
| -VBAT         | TP506 | -VBAT         |
| -VBAT         | TP507 | -VBAT         |
| +VBAT         | TP508 | +VBAT         |
| +VBAT         | TP509 | -VBAT         |

**Test Points**

|  |   |   |
|--|---|---|
| Comments:  | Laboratory:<br><b>Integrated Actuators Laboratory</b>                                 | Variant:<br>Released  |
| Board Name:<br><b>Plume</b>                                      |  | Project Name:<br><b>Smart Footwear for Diabetic Foot Care</b> |
| Sheet Title:<br>Power - Battery Management                       | File Name:<br>Power_Battery_Management.kicad_sch                                      | Designer:<br>Théo Heng  |
| Sheet Path:<br>/Project Architecture/Power - Battery Management/ | Supervisor:<br>Maël Dagon, Paolo Germano  | Date:<br>2025-04-24      Revision:<br>1.0                     |


# [6] Power - Wireless Receiver & USB



**760308101309 Wireless Charging Coil**

| Properties              | Test condition             | Value | Unit | Tol. |
|-------------------------|----------------------------|-------|------|------|
| Inductance              | L = 125 kHz 10 mA          | 10    | μH   | ±10% |
| Rated Current           | I <sub>R</sub> Δt = 40 K   | 1.5   | A    | max. |
| Q-factor                | Q = 125 kHz 10 mA          | 20    |      | typ. |
| Power Capability        | P = V <sub>DC</sub> × 20 V | 20    | W    | typ. |
| DC Resistance           | R <sub>D</sub> @ 20 °C     | 350   | mΩ   | typ. |
| DC Resistance           | R <sub>D</sub> @ 20 °C     | 450   | mΩ   | max. |
| Self Resonant Frequency | f <sub>res</sub>           | 19    | MHz  | typ. |

**Automotive Highly Integrated Wireless Receiver Qi (WPC v2.0) Compliant Power Supply With Dual P-Channel MOSFETs PowerPath.**



**Suggested layout from BQ51013CQWRHLRQ1 evaluation board.**

|  |  |                                 |  |
|--|--|---------------------------------|--|
| <b>Comments:</b><br>EVM User's Guide: BQ51013C-Q1EVM<br>BQ51013C-Q1 Evaluation Module            | <b>Laboratory:</b><br><b>Integrated Actuators Laboratory</b> <b>Variant:</b><br>Released |                                 |  |
| <b>Board Name:</b><br><b>Plume</b> <b>Project Name:</b><br>Smart Footwear for Diabetic Foot Care |  |                                 |  |
| <b>Sheet Title:</b><br>Power - Wireless Receiver & USB   | <b>File Name:</b><br>Power_Wireless_Receiver_USB.kicad                                   | <b>Designer:</b><br>sdfheo Heng | <b>Date:</b><br>2025-04-24             |
| <b>Sheet Path:</b><br>/Project Architecture/Power - Wireless Receiver & USB/                     | <b>Supervisor:</b><br>Maël Dagon, Paolo Germano  | <b>Size:</b><br><b>A3</b>       | <b>Sheet:</b><br><b>6</b> of <b>31</b> |

# [7] Power - Path & Generation

Design Note - Simple Power-Path with PMOS and Schottky diode:


When no USB-C supply nor coil power is present, this simple circuit will by default conduct the battery to the LDO, through the PMOS.

But if USB-C supply or coil power is active, the PMOS will turn off, and the LDO will use the power from the external power supply through the schottky diode.

Design Note - Power Path Schottky diode:

A diode is required to prevent reverse current from flowing to the power source. Selecting the right diode can minimize the leakage current and the forward voltage drop from the power source to the system load. A schottky diode, which has lower forward voltage drop, is recommended. Forward voltage @ If = 200mA: 420mV Reverse current @ Vr = 5V: 70nA


This pull down resistor ensures that the gate is properly pulled to GND when VBUS and VCOIL are not connected.



Active-low reset output, RESET has a push-pull output drive and is capable of directly driving input pins. RESET is low as long as VDD remains below the factory threshold voltage of **3.08V**, and until the delay time ( $t_d = 200\text{ms}$ ) elapses after VDD rises above the threshold voltage.


Current analysis:

**Normal:**  
TPS3839: 150nA  
TPS7A02: 25nA  
PMEG60T20ELR: 65nA  
  
**Shutdown:**  
TPS3839: 150nA  
TPS7A02: 3nA  
PMEG60T20ELR: 60nA



3.3V current draw estimation:  
  
STM32 absolute max: <= 130mA  
Pull-ups: 33uA/Pull-up -> max 1mA  
I2C MUX: <1mA  
I2C Pull-ups: 2.2mA x 2 ~5mA  
Fault LEDs: 2mA/LED -> max 36mA  
RGB LED: (1mA + 3x5mA)x1.25 = 20mA  
IMU: max 1mA  
  
Total Max: 194mA < 200mA

Nanopower IQ, 25-nA, 200-mA, Low-Dropout Voltage Regulator With Fast Transient Response paired with 150-nA, Ultralow Power, Supply Voltage Monitor



Test Points

|   |   |   |
|---|---|---|
| Comments:   | Laboratory:<br><b>Integrated Actuators Laboratory</b> | Variant:<br>Released  |
|   | Board Name:<br><b>Plume</b>                           | Project Name:<br><b>Smart Footwear for Diabetic Foot Care</b> |
| Sheet Title:<br>Power - Path & Generation                       | File Name:<br>Power_Path_Generation.kicad_sch         | Designer:<br>Théo Heng  |
| Sheet Path:<br>/Project Architecture/Power - Path & Generation/ | Supervisor:<br>Maël Dagon, Paolo Germano              | Date:<br>2025-04-24   |
|   |   | Revision:<br>1.0  |
|   | Size:<br><b>A4</b>                                    | Sheet:<br><b>7</b> of <b>31</b>                               |

# [8] Logic - I2C Mux

## A DESIGN NOTE - I2C pull-up resistors:

$R_{MIN} = (VDD_{MAX} - VOL_{MAX}) / IOL$   
 System I/O voltage: VDD = 3.3V +- 5% => 3.47V  
 Low level output voltage (I2C specs): VOL<sub>MAX</sub> = 0.4V  
 Low level output current (I2C specs): IOL = 3mA  
 $R_{MIN} = (3.47V - 0.4V) / 3mA \approx 1k\Omega$

$R_{MAX} \approx (1.18 \times t_{RMAX}) / C_{bMAX}$   
 Standard mode (I2C specs):  $t_{RMAX} = 1000\text{ns}$   
 Fast mode (I2C specs):  $t_{RMAX} = 300\text{ns}$

The maximum bus capacitance for an I2C bus must not exceed 400 pF for fast-mode operation. The bus capacitance can be approximated by adding the capacitance of the STM32,  $C_{STM}$ , the TCA9546A,  $C_{TCA}$ , the capacitance of wires/connections/traces,  $C_{TRACE}$ , and the capacitance of each individual slave (driver) on a given channel  $C_{DRV}$ . If multiple channels will be activated simultaneously, each of the slaves on all channels will contribute to total bus capacitance.



$C_{STM}$  worst case = (not specified in the datasheet, we assume 10pF x 2 = 20pF)  
 $C_{TCA\_IN}$  worst case = 19pF x 2 = 38pF on the SCL, SDA lines  
 $C_{TCA\_OUT}$  worst case = 8pF x 2 = 16pF on the SC3-SC0, SD3-SD0 lines  
 $C_{DRV}$  worst case = (not specified in the datasheet, we assume 8pF x 2 = 16pF)  
 $C_{TRACE}$  = nb of traces x (length of a trace / 25mm) X 2.5pF  
 $C_{TRACE}$  = 2 x (75mm / 25mm) X 2.5pF = 15pF

For one active channel:  
 $C_{bMAX} = C_{STM} + C_{TCA\_IN} + C_{TCA\_OUT} + 9 \times C_{DRV} + C_{TRACE}$   
 $C_{bMAX} = 20pF + 38pF + 16pF + 9 \times 16pF + 15pF \approx 233pF$   
 $\rightarrow R_{MAX} = (1.18 \times 300\text{ns}) / 113pF \approx 1519\Omega$

To be safe, a value of 1k5Ω is chosen. Since both the master I2C lines and fanned channels need pull ups, there will be 2 resistors in parallel, when a channel is activated. Therefore, 3kΩ resistors are chosen, in order to get 1k5Ω.

For two active channel:  
 $C_{bMAX} = C_{STM} + C_{TCA\_IN} + 2x C_{TCA\_OUT} + 18 \times C_{DRV} + 1.66 \times C_{TRACE}$   
 $C_{bMAX} = 20pF + 38pF + 32pF + 18 \times 16pF + 25pF \approx 403pF > 400$

Therefore, only one channel can be enabled at a time, meaning that 9 drivers only can be driven simultaneously.  
 Smaller resistors decrease rise time but increase power consumption. In our case:  $I = 3V3 / 1500 \approx 2.2\text{mA}$




Low Voltage 4-Channel I2C and SMBus Switch with Reset Function

## D Current analysis:

Normal:  
 PCA9546A = 3uA  
 Pull-ups =  $(3.3/1500) \times 2 \times 0.5 = 2.2\text{mA}$   
 Standby:  
 PCA9546A = 1uA  
 Pull-ups = 0uA

MUX\_RESET — TP801 — MUX\_RESET

Test Points



## E Comments:

Laboratory:  
**Integrated Actuators Laboratory**



Variant:  
 Released

Board Name:  
**Plume**



Project Name:  
**Smart Footwear for Diabetic Foot Care**

Sheet Title:  
 Logic - I2C Mux

File Name:  
 Logic\_I2C\_Mux.kicad\_sch

Designer:  
 Théo Heng

Date:  
 2025-04-24

Revision:  
 1.0



Sheet Path:  
 /Project Architecture/Logic - I2C Mux/

Supervisor:  
 Maël Dagon, Paolo Germano

Size:  
**A4**

Sheet:  
**8** of **31**


# [9] User - RGB LED Indicator



**Current analysis:**

Normal:  
REG71050: 65µA  
SN74AHCT1G125: 10µA  
IN-PI22TAT5R5G5B: 1mA

Shutdown:  
REG71050: 10nA  
SN74AHCT1G125: OFF  
IN-PI22TAT5R5G5B: OFF




**Comments:**  
This LED was chosen as it is relatively small compared to other well-known alternatives (e.g. WS2812B). The saved space is used to provide proper 5V supply and logic. This LED might work with 3V3 but this approach is safer.

|   |                                  |                                       |
|---|----------------------------------|---------------------------------------|
| Laboratory:                                     | <b>EPFL</b>                      | Variant:                              |
|   | Plume                            | Released                              |
| Board Name:                                     | <b>Plume</b>                     | Project Name:                         |
|   |                                  | Smart Footwear for Diabetic Foot Care |
| Sheet Title:                                    | File Name:                       | Date:                                 |
| User - RGB LED Indicator                        | User_RGB_LED_Indicator.kicad_sch | Théo Heng<br>2025-04-24               |
| Sheet Path:                                     | Designer:                        | Revision:                             |
| /Project Architecture/User - RGB LED Indicator/ | Maël Dagon, Paolo Germano        | 1.0                                   |
| Supervisor:                                     | Size:                            | Sheet:                                |
|   | <b>A4</b>                        | <b>9</b> of <b>31</b>                 |

# [10] Sensing - IMU

IMU\_INT1 → TP1001 IMU\_INT1  
IMU\_INT2 → TP1002 IMU\_INT2

Test Points



6-axis, smart, low power Inertial Measurement Unit for high-performance applications



LGA-14 2.5x3x0.83 mm




Current analysis BMI270:  
Normal: A+G Performance mode: 970uA  
Wake-up mode: 5uA  
Sleep: 3uA

DESIGN NOTE - I<sup>2</sup>C address:

The default I<sup>2</sup>C address of the device is 0b1101000 (0x68). It is used if the SDO pin is pulled to 'GND'. The alternative address 0b1101001 (0x69) is selected by pulling the SDO pin to 'VDDIO'.

|                                      |   |           |                 |
|--------------------------------------|---|-----------|-----------------|
| Comments:                            | Laboratory:                            |           | Variant:        |
|                                      |   |           | Released        |
| Board Name:                          | Project Name:   |           |                 |
| Plume                                |  Smart Footwear for Diabetic Foot Care |           |                 |
| Sheet Title:                         | File Name:  | Designer: | Date: Revision: |
| Sensing - IMU                        | Sensing_IMU.kicad_sch   | Théo Heng | 2025-04-24 1.0  |
| Sheet Path:                          | Supervisor:   | Size:     | Sheet:          |
| /Project Architecture/Sensing - IMU/ | Maël Dagon, Paolo Germano   | A4        | 10 of 31        |

# [11] Motor Control - Drivers x18



**30 positions with screw: 44.48mm wide  
25 positions with screw: 39.48mm wide  
15 positions with screw: 29.48mm wide**

**DEBUG NOTE:**  
The specific Samtec connectors used on this PCB are not usually in stock on Mouser/Digikey so they are directly ordered from Samtech's website.

**LAYOUT NOTE:**  
Use the "Replicate Layout" KiCad Plugin to layout copy-paste the layout of each driver unit.



**DESIGN NOTE:**  
This first iteration implements 18 drivers module, with 11 of them used at the front of the foot, and 7 at the rear.  
  
Each module is similar, with the DRV8214RTER IC, decoupling and bulk capacitors, test points, fault LED, and current sense resistor.  
  
For the layout, every module are routed in the exact same way.

nSLEEP\_FRONT → TP1101 nSLEEP\_FRONT  
nSLEEP\_REAR → TP1102 nSLEEP\_REAR  
SD0 → TP1103 SD0  
SC0 → TP1104 SC0  
SD1 → TP1105 SD1  
SC1 → TP1106 SC1

## Test Points



|   |   |                                      |
|---|---|--------------------------------------|
| Comments:   | Laboratory:<br><b>Integrated Actuators Laboratory</b>         | Variant:<br>Released                 |
| Board Name:<br><b>Plume</b>                                   | Project Name:<br><b>Smart Footwear for Diabetic Foot Care</b> |                                      |
| Sheet Title:<br>Motor Control - Drivers                       | File Name:<br>Motor_Control_Drivers.kicad_sch                 | Designer:<br>Théo Heng               |
| Sheet Path:<br>/Project Architecture/Motor Control - Drivers/ | Supervisor:<br>Maël Dagon, Paolo Germano                      | Date:<br>2025-04-24 Revision:<br>1.0 |

# [12] Motor Control - Driver 12





|   |  |                         |
|---|--|-------------------------|
| Comments:   | Laboratory:                                  | Variant:                |
|   | <b>EPFL</b><br><b>Plume</b>                  | Released                |
| Board Name:   | Project Name:                                |                         |
| <b>Plume</b>  | <b>Smart Footwear for Diabetic Foot Care</b> |                         |
| Sheet Title:  | File Name:                                   | Date:                   |
| Motor Control - Driver  | Motor_Control_Driver.kicad_sch               | Théo Heng<br>2025-04-24 |
| Sheet Path:   | Designer:                                    | Revision:               |
| /Project Architecture/Motor Control - Drivers/Motor Control - Driver 0/ | Maël Dagon, Paolo Germano                    | 1.0                     |
| Supervisor:   | Size:  | Sheet:                  |
|   | <b>A4</b>                                    | <b>12 of 31</b>         |

# [13] Motor Control - Driver 13





|   |                                |                                       |
|---|--------------------------------|---------------------------------------|
| Comments:   | Laboratory:                    | Variant:                              |
|   | <b>EPFL</b>                    | Released                              |
| Board Name:   | Plume                          | Project Name:                         |
|   | <i>Plume</i>                   | Smart Footwear for Diabetic Foot Care |
| Sheet Title:  | File Name:                     | Date:                                 |
| Motor Control - Driver  | Motor_Control_Driver.kicad_sch | Théo Heng                             |
| Sheet Path:   | Designer:                      | Revision:                             |
| /Project Architecture/Motor Control - Drivers/Motor Control - Driver 1/ | Supervisor:                    | 2025-04-24                            |
|   | Maël Dagon, Paolo Germano      | 1.0                                   |
|   | Size:                          | Sheet:                                |
|   | <b>A4</b>                      | <b>13</b> of <b>31</b>                |

# [14] Motor Control - Driver 14





|  |  |   |
|--|--|---|
| Comments:  | Laboratory:                                  | Variant:  |
|  | <b>EPFL</b>                                  | Released  |
| Board Name:<br><b>Plume</b>  | <b>Plume</b>                                 | Project Name:<br><b>Smart Footwear for Diabetic Foot Care</b> |
| Sheet Title:<br>Motor Control - Driver   | File Name:<br>Motor_Control_Driver.kicad_sch | Designer:<br>Théo Heng  |
| Sheet Path:<br>/Project Architecture/Motor Control - Drivers/Motor Control - Driver 2/ | Supervisor:<br>Maël Dagon, Paolo Germano     | Date:<br>2025-04-24   |
|  | Size:<br><b>A4</b>                           | Revision:<br>1.0  |

# [15] Motor Control - Driver 15



|  |  |   |
|--|--|---|
| Comments:  | Laboratory:                                  | Variant:  |
|  | <b>EPFL</b>                                  | Released  |
| Board Name:<br><b>Plume</b>  | <b>Plume</b>                                 | Project Name:<br><b>Smart Footwear for Diabetic Foot Care</b> |
| Sheet Title:<br>Motor Control - Driver   | File Name:<br>Motor_Control_Driver.kicad_sch | Designer:<br>Théo Heng  |
| Sheet Path:<br>/Project Architecture/Motor Control - Drivers/Motor Control - Driver 3/ | Supervisor:<br>Maël Dagon, Paolo Germano     | Date: 2025-04-24 Revision: 1.0                                |
|  | Size:<br><b>A4</b>                           | Sheet:<br><b>15</b> of <b>31</b>                              |

# [16] Motor Control - Driver 16



|  |  |   |
|--|--|---|
| Comments:  | Laboratory:                                  | Variant:  |
|  | <b>EPFL</b>                                  | Released  |
| Board Name:<br><b>Plume</b>  | <b>Plume</b>                                 | Project Name:<br><b>Smart Footwear for Diabetic Foot Care</b> |
| Sheet Title:<br>Motor Control - Driver   | File Name:<br>Motor_Control_Driver.kicad_sch | Designer:<br>Théo Heng  |
| Sheet Path:<br>/Project Architecture/Motor Control - Drivers/Motor Control - Driver 4/ | Supervisor:<br>Maël Dagon, Paolo Germano     | Date:<br>2025-04-24   |
|  | Size:<br><b>A4</b>                           | Revision:<br>1.0  |

# [17] Motor Control - Driver 17

**A DESIGN NOTE - Bulk Capacitance:**

Having appropriate local bulk capacitance is an important factor in motor drive system design. Having more bulk capacitance is generally beneficial, while the disadvantages are increased cost and physical size.

The amount of local capacitance needed depends on a variety of factors, including:

- The highest current required by the motor system
- The capacitance of the power supply and ability to source current
- The amount of parasitic inductance between the power supply and motor system
- The acceptable voltage ripple
- The type of motor used (brushed DC, brushless DC, stepper)
- The motor braking method

The inductance between the power supply and motor drive system limits how the rate current can change from the power supply. If the local bulk capacitance is too small, the system responds to excessive current demands or dumps from the motor with a change in voltage. When adequate bulk capacitance is used, the motor voltage remains stable and high current can be quickly supplied.

The data sheet generally provides a recommended value, but system-level testing is required to determine the appropriate sized bulk capacitor.

The voltage rating for bulk capacitors should be higher than the operating voltage, to provide margin for cases when the motor transfers energy to the supply.

**B DESIGN NOTE - IPROPI Resistor:**

The current regulation threshold (I<sub>TRIP</sub>) is set through a combination of the V<sub>REF</sub> voltage and I<sub>PROPI</sub> output resistor (R<sub>PROPI</sub>). This is done by comparing the voltage drop across the external R<sub>PROPI</sub> resistor to V<sub>REF</sub> with an internal comparator

$$I_{TRIP} (A) = V_{REF} (V) / (R_{PROPI} (\Omega) \times A_{PROPI} (\mu A/A))$$

$$I_{TRIP} (A) = 0.5V / (680\Omega \times 5560\mu A/A) \approx 132mA$$



**C DESIGN NOTE - FAULT LED:**

APTD1608LQBC/D has a typical luminous intensity of 20mcd (millicandela) @2mA  
Typ. V<sub>f</sub> = 2.65V @2mA,

$$Typ. I<sub>f</sub> = (3.3V - 2.65V)/330hm \approx 2 mA$$

$$P = 2.65V \times 2mA \approx 5 mW < 120 mW (max rating)$$

**DRV8214 2-A Brushed DC Motor Driver with Ripple Counting, Stall Detection, and Speed Regulation**



WQFN-16 3x3 mm

Current analysis DRV8214:



Normal mode motor stopped: 2.8 mA  
Sleep: 100nA

|        |        |        |
|--------|--------|--------|
| nFAULT | TP1701 | nFAULT |
| RC_OUT | TP1702 | RC_OUT |
| OUT1   | TP1703 | OUT1   |
| OUT2   | TP1704 | OUT2   |
| A0     | TP1705 | A0     |
| A1     | TP1706 | A1     |

## Test Points



|  |  |   |
|--|--|---|
| Comments:  | Laboratory:                                  | Variant:  |
|  | <b>EPFL</b>                                  | Released  |
| Board Name:<br><b>Plume</b>  | <b>Plume</b>                                 | Project Name:<br><b>Smart Footwear for Diabetic Foot Care</b> |
| Sheet Title:<br>Motor Control - Driver   | File Name:<br>Motor_Control_Driver.kicad_sch | Designer:<br>Théo Heng  |
| Sheet Path:<br>/Project Architecture/Motor Control - Drivers/Motor Control - Driver 5/ | Supervisor:<br>Maël Dagon, Paolo Germano     | Date:<br>2025-04-24   |
|  | Size:<br><b>A4</b>                           | Revision:<br>1.0  |
|  | Sheet:<br><b>17</b>                          | of <b>31</b>  |

# [18] Motor Control - Driver 18





|  |  |   |
|--|--|---|
| Comments:  | Laboratory:                                  | Variant:  |
|  | <b>EPFL</b>                                  | Released  |
| Board Name:<br><b>Plume</b>  | <b>Plume</b>                                 | Project Name:<br><b>Smart Footwear for Diabetic Foot Care</b> |
| Sheet Title:<br>Motor Control - Driver   | File Name:<br>Motor_Control_Driver.kicad_sch | Designer:<br>Théo Heng  |
| Sheet Path:<br>/Project Architecture/Motor Control - Drivers/Motor Control - Driver 6/ | Supervisor:<br>Maël Dagon, Paolo Germano     | Date:<br>2025-04-24   |
|  | Size:<br><b>A4</b>                           | Revision:<br>1.0  |

# [19] Motor Control - Driver 19



|  |  |   |
|--|--|---|
| Comments:  | Laboratory:                                  | Variant:  |
|  | <b>EPFL</b>                                  | Released  |
| Board Name:<br><b>Plume</b>  | <b>Plume</b>                                 | Project Name:<br><b>Smart Footwear for Diabetic Foot Care</b> |
| Sheet Title:<br>Motor Control - Driver   | File Name:<br>Motor_Control_Driver.kicad_sch | Designer:<br>Théo Heng  |
| Sheet Path:<br>/Project Architecture/Motor Control - Drivers/Motor Control - Driver 7/ | Supervisor:<br>Maël Dagon, Paolo Germano     | Date:<br>2025-04-24   |
|  | Size:<br><b>A4</b>                           | Revision:<br>1.0  |

# [20] Motor Control - Driver 20



|  |   |                        |
|--|---|------------------------|
| Comments:  | Laboratory:   | Variant:               |
|  | <b>EPFL</b>   | Released               |
| Board Name:<br><b>Plume</b>  | Project Name:<br><b>Smart Footwear for Diabetic Foot Care</b> |                        |
| Sheet Title:<br>Motor Control - Driver   | File Name:<br>Motor_Control_Driver.kicad_sch                  | Designer:<br>Théo Heng |
| Sheet Path:<br>/Project Architecture/Motor Control - Drivers/Motor Control - Driver 8/ | Supervisor:<br>Maël Dagon, Paolo Germano                      | Date:<br>2025-04-24    |
|  | Size:<br><b>A4</b>  | Revision:<br>1.0       |

[21] Motor Control - Driver 21

## DESIGN NOTE - Bulk Capacitance:

Having appropriate local bulk capacitance is an important factor in motor drive system design. Having more bulk capacitance is generally beneficial, while the disadvantages are increased cost and physical size.

The amount of local capacitance needed depends on a variety of factors, including:

- The highest current required by the motor system
  - The capacitance of the power supply and ability to source current
  - The amount of parasitic inductance between the power supply and motor system
  - The acceptable voltage ripple
  - The type of motor used (brushed DC, brushless DC, stepper)
  - The motor braking method

The inductance between the power supply and motor drive system limits how the rate current can change from the power supply. If the local bulk capacitance is too small, the system responds to excessive current demands or dumps from the motor with a change in voltage. When adequate bulk capacitance is used, the motor voltage remains stable and high current can be quickly supplied.

The data sheet generally provides a recommended value, but system-level testing is required to determine the appropriate sized bulk capacitor.

The voltage rating for bulk capacitors should be higher than the operating voltage, to provide margin for cases when the motor transfers energy to the supply.

DESIGN NOTE - IPROPI Resisto

The current regulation threshold (ITRIP) is set through a combination of the VREF voltage and IPROPI output resistor (RIPROPI). This is done by comparing the voltage drop across the external RIPROPI resistor to VVREF with an internal comparator



$$ITRIP(A) = VVREF(V) / (RIPROPI(\Omega) \times AIPROPI(\mu A/A))$$

DESIGN NOTE - FAULT LED

APTD1608LQBC/D has a typical luminous intensity of 20mcd (millicandela) @2mA  
Typ. Vf = 2.65V @2mA.

$$\text{Typ. } I_f = (3.3V - 2.65V)/330\text{hm} \approx 2 \text{ mA}$$

$$P = 2.65V \times 2\text{mA} \approx 5 \text{ mW} < 120 \text{ mW (max rating)}$$



WQFN-16 3x3 mm

Current analysis DRV8214:

Normal mode motor stopped: 2.8 mA



Sleep: 100nA

## **DRV8214 2-A Brushed DC Motor Driver with Ripple Counting, Stall Detection, and Speed Regulation**

nFAULT —○— TP2101 nFAULT  
 RC\_OUT —○— TP2102 RC\_OUT  
 OUT1 —○— TP2103 OUT1  
 OUT2 —○— TP2104 OUT2  
 A0 —○— TP2105 A0  
 A1 —○— TP2106 A1



|  |  |  |
|--|--|--|
| Comments:  | Laboratory:<br><br> | Variant:<br>Released                                   |
|  | Board Name:<br><br> | Project Name:<br>Smart Footwear for Diabetic Foot Care |
| Sheet Title:<br>Motor Control - Driver   | File Name:<br>Motor_Control_Driver.kicad_sch   | Designer:<br>Théo Heng                                 |
|  |  | Date: 2025-04-24 Revision: 1.0                         |
| Sheet Path:<br>/Project Architecture/Motor Control - Drivers/Motor Control - Driver 9/ | Supervisor:<br>Maël Dagon, Paolo Germano   | Size: A4 Sheet: 21 of 31                               |

# [22] Motor Control - Driver 22



|  |                                |                                       |
|--|--------------------------------|---------------------------------------|
| Comments:  | Laboratory:                    | Variant:                              |
|  | <b>EPFL</b>                    | Released                              |
| Board Name:  | <b>Plume</b>                   | Project Name:                         |
| Motor Control - Driver   |                                | Smart Footwear for Diabetic Foot Care |
| Sheet Title:   | File Name:                     | Date:                                 |
| Motor Control - Driver   | Motor_Control_Driver.kicad_sch | Théo Heng                             |
| Sheet Path:  | Designer:                      | Revision:                             |
| /Project Architecture/Motor Control - Drivers/Motor Control - Driver 10/ | Maël Dagon, Paolo Germano      | 2025-04-24                            |
|  | Supervisor:                    | Size:                                 |
|  | Maël Dagon, Paolo Germano      | A4                                    |
|  | Sheet:                         | 22 of 31                              |

# [23] Motor Control - Driver 23



|   |  |   |
|---|--|---|
| Comments:   | Laboratory:                                  | Variant:  |
|   | <b>EPFL</b>                                  | Released  |
| Board Name:<br><b>Plume</b>   | <b>Plume</b>                                 | Project Name:<br><b>Smart Footwear for Diabetic Foot Care</b> |
| Sheet Title:<br>Motor Control - Driver  | File Name:<br>Motor_Control_Driver.kicad_sch | Designer:<br>Théo Heng  |
| Sheet Path:<br>/Project Architecture/Motor Control - Drivers/Motor Control - Driver 11/ | Supervisor:<br>Maël Dagon, Paolo Germano     | Date:<br>2025-04-24   |
|   | Size:<br><b>A4</b>                           | Revision:<br>1.0  |

# [24] Motor Control - Driver 24

**DESIGN NOTE - Bulk Capacitance:**

Having appropriate local bulk capacitance is an important factor in motor drive system design. Having more bulk capacitance is generally beneficial, while the disadvantages are increased cost and physical size.

The amount of local capacitance needed depends on a variety of factors, including:

- The highest current required by the motor system
- The capacitance of the power supply and ability to source current
- The amount of parasitic inductance between the power supply and motor system
- The acceptable voltage ripple
- The type of motor used (brushed DC, brushless DC, stepper)
- The motor braking method

The inductance between the power supply and motor drive system limits how the rate current can change from the power supply. If the local bulk capacitance is too small, the system responds to excessive current demands or dumps from the motor with a change in voltage. When adequate bulk capacitance is used, the motor voltage remains stable and high current can be quickly supplied.

The data sheet generally provides a recommended value, but system-level testing is required to determine the appropriate sized bulk capacitor.

The voltage rating for bulk capacitors should be higher than the operating voltage, to provide margin for cases when the motor transfers energy to the supply.

**DESIGN NOTE - IPROPI Resistor:**

The current regulation threshold (I<sub>TRIP</sub>) is set through a combination of the V<sub>REF</sub> voltage and I<sub>PROPI</sub> output resistor (R<sub>PROPI</sub>). This is done by comparing the voltage drop across the external R<sub>PROPI</sub> resistor to V<sub>REF</sub> with an internal comparator

$$I_{TRIP} (A) = V_{REF} (V) / (R_{PROPI} (Ω) × A_{PROPI} (\mu A/A))$$



$$I_{TRIP} (A) = 0.5V / (680Ω × 5560\mu A/A) ≈ 132mA$$

**DESIGN NOTE - FAULT LED:**

APTD1608LQBC/D has a typical luminous intensity of 20mcd (millilumen) @2mA  
Typ. V<sub>f</sub> = 2.65V @2mA,

$$Typ. I<sub>f</sub> = (3.3V - 2.65V)/330Ω ≈ 2 mA$$

$$P = 2.65V × 2mA ≈ 5 mW < 120 mW (max rating)$$


**WQFN-16 3x3 mm**
**Current analysis DRV8214:**

Normal mode motor stopped: 2.8 mA  
Sleep: 100nA



**DRV8214 2-A Brushed DC Motor Driver with Ripple Counting, Stall Detection, and Speed Regulation**

|        |        |        |
|--------|--------|--------|
| nFAULT | TP2401 | nFAULT |
| RC_OUT | TP2402 | RC_OUT |
| OUT1   | TP2403 | OUT1   |
| OUT2   | TP2404 | OUT2   |
| A0     | TP2405 | A0     |
| A1     | TP2406 | A1     |

**Test Points**



|   |  |   |
|---|--|---|
| Comments:   | Laboratory:                                  | Variant:  |
|   | <b>EPFL</b>                                  | Released  |
| Board Name:<br><b>Plume</b>   | <b>Plume</b>                                 | Project Name:<br><b>Smart Footwear for Diabetic Foot Care</b> |
| Sheet Title:<br>Motor Control - Driver  | File Name:<br>Motor_Control_Driver.kicad_sch | Designer:<br>Théo Heng  |
| Sheet Path:<br>/Project Architecture/Motor Control - Drivers/Motor Control - Driver 12/ | Supervisor:<br>Maël Dagon, Paolo Germano     | Date:<br>2025-04-24   |
|   | Size:<br><b>A4</b>                           | Revision:<br>1.0  |

# [25] Motor Control - Driver 25



|  |  |                         |
|--|--|-------------------------|
| Comments:  | Laboratory:                                  | Variant:                |
|  | <b>EPFL</b><br><b>Plume</b>                  | Released                |
| Board Name:  | Project Name:                                |                         |
| <b>Plume</b>   | <b>Smart Footwear for Diabetic Foot Care</b> |                         |
| Sheet Title:   | File Name:                                   | Date:                   |
| Motor Control - Driver   | Motor_Control_Driver.kicad_sch               | Théo Heng<br>2025-04-24 |
| Sheet Path:  | Designer:                                    | Revision:               |
| /Project Architecture/Motor Control - Drivers/Motor Control - Driver 13/ | Maël Dagon, Paolo Germano                    | 1.0                     |
| Supervisor:  | Size:  | Sheet:                  |
|  | <b>A4</b>                                    | <b>25</b> of <b>31</b>  |

# [26] Motor Control - Driver 26



|   |  |   |
|---|--|---|
| Comments:   | Laboratory:                                  | Variant:  |
|   | <b>EPFL</b>                                  | Released  |
| Board Name:<br><b>Plume</b>   | <b>Plume</b>                                 | Project Name:<br><b>Smart Footwear for Diabetic Foot Care</b> |
| Sheet Title:<br>Motor Control - Driver  | File Name:<br>Motor_Control_Driver.kicad_sch | Designer:<br>Théo Heng  |
| Sheet Path:<br>/Project Architecture/Motor Control - Drivers/Motor Control - Driver 14/ | Supervisor:<br>Maël Dagon, Paolo Germano     | Date:<br>2025-04-24   |
|   | Size:<br><b>A4</b>                           | Revision:<br>1.0  |

# [27] Motor Control - Driver 27

**DESIGN NOTE - Bulk Capacitance:**

Having appropriate local bulk capacitance is an important factor in motor drive system design. Having more bulk capacitance is generally beneficial, while the disadvantages are increased cost and physical size.

The amount of local capacitance needed depends on a variety of factors, including:

- The highest current required by the motor system
- The capacitance of the power supply and ability to source current
- The amount of parasitic inductance between the power supply and motor system
- The acceptable voltage ripple
- The type of motor used (brushed DC, brushless DC, stepper)
- The motor braking method

The inductance between the power supply and motor drive system limits how the rate current can change from the power supply. If the local bulk capacitance is too small, the system responds to excessive current demands or dumps from the motor with a change in voltage. When adequate bulk capacitance is used, the motor voltage remains stable and high current can be quickly supplied.

The data sheet generally provides a recommended value, but system-level testing is required to determine the appropriate sized bulk capacitor.

The voltage rating for bulk capacitors should be higher than the operating voltage, to provide margin for cases when the motor transfers energy to the supply.

**DESIGN NOTE - IPROPI Resistor:**

The current regulation threshold (I<sub>TRIP</sub>) is set through a combination of the V<sub>REF</sub> voltage and I<sub>PROPI</sub> output resistor (R<sub>PROPI</sub>). This is done by comparing the voltage drop across the external R<sub>PROPI</sub> resistor to V<sub>REF</sub> with an internal comparator

$$I_{TRIP} (A) = V_{REF} (V) / (R_{PROPI} (\Omega) \times A_{PROPI} (\mu A/A))$$



$$I_{TRIP} (A) = 0.5V / (680\Omega \times 5560\mu A/A) \approx 132mA$$

**DESIGN NOTE - FAULT LED:**

APTD1608LQBC/D has a typical luminous intensity of 20mcd (millicandela) @2mA  
Typ. V<sub>f</sub> = 2.65V @2mA,

$$Typ. I<sub>f</sub> = (3.3V - 2.65V)/330hm \approx 2 mA$$

$$P = 2.65V \times 2mA \approx 5 mW < 120 mW (max rating)$$


**WQFN-16 3x3 mm**
**Current analysis DRV8214:**

Normal mode motor stopped: 2.8 mA  
Sleep: 100nA



**DRV8214 2-A Brushed DC Motor Driver with Ripple Counting, Stall Detection, and Speed Regulation**

|        |        |        |
|--------|--------|--------|
| nFAULT | TP2701 | nFAULT |
| RC_OUT | TP2702 | RC_OUT |
| OUT1   | TP2703 | OUT1   |
| OUT2   | TP2704 | OUT2   |
| A0     | TP2705 | A0     |
| A1     | TP2706 | A1     |

**Test Points**

|   |  |   |
|---|--|---|
| Comments:   | Laboratory:                                  | Variant:  |
|   | <b>EPFL</b>                                  | Released  |
| Board Name:<br><b>Plume</b>   | <b>Plume</b>                                 | Project Name:<br><b>Smart Footwear for Diabetic Foot Care</b> |
| Sheet Title:<br>Motor Control - Driver  | File Name:<br>Motor_Control_Driver.kicad_sch | Designer:<br>Théo Heng  |
| Sheet Path:<br>/Project Architecture/Motor Control - Drivers/Motor Control - Driver 15/ | Supervisor:<br>Maël Dagon, Paolo Germano     | Date:<br>2025-04-24   |
|   | Size:<br><b>A4</b>                           | Revision:<br>1.0  |

# [28] Motor Control - Driver 28



|   |   |                        |
|---|---|------------------------|
| Comments:   | Laboratory:   | Variant:               |
|   | <b>EPFL</b>   | Released               |
| Board Name:<br><b>Plume</b>   | Project Name:<br><b>Smart Footwear for Diabetic Foot Care</b> |                        |
| Sheet Title:<br>Motor Control - Driver  | File Name:<br>Motor_Control_Driver.kicad_sch                  | Designer:<br>Théo Heng |
| Sheet Path:<br>/Project Architecture/Motor Control - Drivers/Motor Control - Driver 16/ | Supervisor:<br>Maël Dagon, Paolo Germano                      | Date:<br>2025-04-24    |
|   | Size:<br><b>A4</b>  | Revision:<br>1.0       |

# [29] Motor Control - Driver 29

**DESIGN NOTE - Bulk Capacitance:**  
Having appropriate local bulk capacitance is an important factor in motor drive system design. Having more bulk capacitance is generally beneficial, while the disadvantages are increased cost and physical size.

The amount of local capacitance needed depends on a variety of factors, including:  
 • The highest current required by the motor system  
 • The capacitance of the power supply and ability to source current  
 • The amount of parasitic inductance between the power supply and motor system  
 • The acceptable voltage ripple  
 • The type of motor used (brushed DC, brushless DC, stepper)  
 • The motor braking method

The inductance between the power supply and motor drive system limits how the rate current can change from the power supply. If the local bulk capacitance is too small, the system responds to excessive current demands or dumps from the motor with a change in voltage. When adequate bulk capacitance is used, the motor voltage remains stable and high current can be quickly supplied.

The data sheet generally provides a recommended value, but system-level testing is required to determine the appropriate sized bulk capacitor.



The voltage rating for bulk capacitors should be higher than the operating voltage, to provide margin for cases when the motor transfers energy to the supply.

**DESIGN NOTE - I<sub>PROPI</sub> Resistor:**  
The current regulation threshold (I<sub>TRIP</sub>) is set through a combination of the V<sub>REF</sub> voltage and I<sub>PROPI</sub> output resistor (R<sub>I<sub>PROPI</sub></sub>). This is done by comparing the voltage drop across the external R<sub>I<sub>PROPI</sub></sub> resistor to V<sub>REF</sub> with an internal comparator

$$I_{TRIP} (A) = V_{REF} (V) / (R_{I\text{PROPI}} (\Omega) \times A_{I\text{PROPI}} (\mu A/A))$$

$$I_{TRIP} (A) = 0.5V / (680\Omega \times 5560\mu A/A) \approx 132mA$$

**DESIGN NOTE - FAULT LED:**  
APTD1608LQBC/D has a typical luminous intensity of 20mcd (millilumen) @2mA  
Typ.  $V_f = 2.65V$  @2mA,  
 $I_f = (3.3V - 2.65V)/330\Omega \approx 2 mA$   
 $P = 2.65V \times 2mA \approx 5 mW < 120 mW$  (max rating)



WQFN-16 3x3 mm

Current analysis DRV8214:

Normal mode motor stopped: 2.8 mA  
Sleep: 100nA

DRV8214 2-A Brushed DC Motor Driver with Ripple Counting, Stall Detection, and Speed Regulation

|        |        |        |
|--------|--------|--------|
| nFAULT | TP2901 | nFAULT |
| RC_OUT | TP2902 | RC_OUT |
| OUT1   | TP2903 | OUT1   |
| OUT2   | TP2904 | OUT2   |
| A0     | TP2905 | A0     |
| A1     | TP2906 | A1     |

## Test Points

|   |  |   |
|---|--|---|
| Comments:   | Laboratory:                                  | Variant:  |
|   | <b>EPFL</b>                                  | Released  |
| Board Name:<br><b>Plume</b>   | <b>Plume</b>                                 | Project Name:<br><b>Smart Footwear for Diabetic Foot Care</b> |
| Sheet Title:<br>Motor Control - Driver  | File Name:<br>Motor_Control_Driver.kicad_sch | Designer:<br>Théo Heng  |
| Sheet Path:<br>/Project Architecture/Motor Control - Drivers/Motor Control - Driver 17/ | Supervisor:<br>Maël Dagon, Paolo Germano     | Date:<br>2025-04-24   |
|   | Size:<br><b>A4</b>                           | Revision:<br>1.0  |


# [30] Misc - Mounting Holes

## DEBUG NOTE:

HOLE 1 is connected to GND and can be very useful if used as a GND reference for multimeter debug.

## DESIGN NOTE:

Fiducials are omitted due to space constraints and are in any case unnecessary in the foreseeable future because manual assembly is required since the DRV8214 ICs, among others, are unavailable on JLCPCB parts store.



Mechanical Mounting holes

- FID3001 Fiducial
- FID3002 Fiducial
- FID3005 Fiducial
- FID3006 Fiducial

Fiducials Front

- H3005 MountingHole
- H3006 MountingHole
- H3007 MountingHole
- H3008 MountingHole

Assembly Frame Mounting Holes

## DESIGN NOTE:

In general, do not ground mounting holes that are connected to chassis if there is a risk that the chassis sends current into the PCB and damages it or if the PCB works with high voltage and there is a risk for the user which is in contact with the chassis.

For a PCB in a metallic enclosure, the enclosure can act like a big antenna. The cavities (space between the ground plane on the PCB and the bottom metallic plane of the enclosure) cause capacitive coupling. This is because we have two metallic planes separated by a dielectric, at a separated potential, causing common mode noise, and result in the noise being amplified by the metal enclosure. Ideally, in this case, the case should be grounded. But the grounding will cause multiple return paths if multiple mounting holes are connected to GND. In that case, some current will flow in the enclosure and can give a little shock when the user touches it. This can be fine for small battery-operated devices.

The best option is therefore to ground only one mounting hole.

- FID3003 Fiducial
- FID3004 Fiducial
- FID3007 Fiducial
- FID3008 Fiducial

Fiducials Rear

|  |   |   |
|--|---|---|
| Comments:<br>Inspired by Amulet controller Schematics<br>by Vincent Nguyen<br><br>How to connect mounting holes.<br>PCB Mounting Holes | Laboratory:<br><br><br><br> | Variant:<br>Released                      |
| Board Name:<br><br><b>Plume</b>  | Project Name:<br><br><b>Smart Footwear for Diabetic Foot Care</b>   |   |
| Sheet Title:<br><br>Misc - Mounting Holes  | File Name:<br><br>Misc_Mounting_Holes.kicad_sch   | Designer:<br><br>Théo Heng                |
| Sheet Path:<br><br>/Project Architecture/Misc - Mounting Holes & Fiducials/  | Supervisor:<br><br>Maël Dagon, Paolo Germano  | Date:<br>2025-04-24      Revision:<br>1.0 |
|  | Size:<br><br><b>A4</b>  | Sheet:<br><br><b>30</b> of <b>31</b>      |

## [31] Revision History

DD.MM.YYYY - xxx Revision  
Variant: V1.0 -> V1.1

Ordered 24/04/2025: 5 PCB for 232.48\$  
Components for 4 PCB on Mouser: 300 CHF

- - Changed Silkscreen for motor testpoints from D10 to D8
  - ✗ - GND test point next to battery connector isn't actually a GND test point -> changed to BAT-
  - ✗ - Make the slot for the battery cables wider, move the connector further inside the PCB.
  - ✗ - Move R & C of holes further in the bottom
  - ✗ - Power LED is not ON when VBUS is present
  - ✗ - Paste mask of the STM32 MCU is bad apparently
  - ✗ - Need to short +VBAT and DOUT/COUT to enable the device
  - ✗ - Switch 1 and switch2 needs to be routed to different GPIO to have them as interrupt

DD.MM.YYYY - xxx Revision  
Variant: xxx

DD.MM.YYYY - xxx Revision  
Variant: xxx

DD.MM.YYYY - xxx Revision  
Variant: xxx

- Done
- ✗ Not done yet

|                                   |   |   |   |                                  |
|-----------------------------------|---|---|---|----------------------------------|
| Comments:                         | Laboratory:<br><b>Integrated Actuators Laboratory</b> | <b>EPFL</b>   | Variant:<br>Released  |                                  |
|                                   | Board Name:<br><b>Plume</b>                           |  | Project Name:<br><b>Smart Footwear for Diabetic Foot Care</b> |                                  |
| Sheet Title:<br>Revision History  | File Name:<br>Revision History.kicad_sch              | Designer:<br>Théo Heng  | Date:<br>2025-04-24   | Revision:<br>1.0                 |
| Sheet Path:<br>/Revision History/ | Superviser:<br>Maël Dagon, Paolo Germano              |   | Size:<br><b>A4</b>  | Sheet:<br><b>31</b> of <b>31</b> |